

Instructions for cleaning and maintenance

The appliances (splints) fabricated out of thermoforming materials should be cleaned and maintained as follows:

- Best results are achieved with **Oxydens** cleansing tablets for dental splints (280 030, Oxydens Clean-set, 280 032, 32 cleansing tablets).
Further cleaning agents: Soap, curd soap, liquid soap and dish-washers.
Do not use any strongly perfumed soaps.
Not suited are: tooth-paste (abrasives), mouth-wash (can cause discolouring) and water hotter than 50 °C (deformation).

After use:

- Well wash with water.
- With no Oxydens at hand thoroughly clean the inner and exterior side of the splint with a tooth brush and soap.
- Again well wash with water.
- Shake off the water or dry with a towel.
- Never blow-dry - deformation!
- Very important: Allow the splint to completely dry!
Keep at a dry place, at best in a box like the Erkobox (215 030) or Splintbox (214 020) that has aeration holes.
- Again wash with water before using it.

Generation of bad smells

If after some time the splint has adopted a strong smell, apply Oxydens or put the splint for one hour in a non-perfumed, concentrated soap solution, afterwards thoroughly wash with water. Such a soap sud removes most of the smell generating bacteria.

Discolouring

Soft thermoforming materials have the tendency to discolour. This intake of colour pigments can be reduced or avoided by a careful maintenance but it cannot be reversed. Mouth-washes and amalgam fillings can also cause discolouring.

Disinfection

With the exception of Erkoloc and Erkoloc-pro all thermoforming materials can be disinfected with disinfection alcohol and other commercial liquids. Erkoloc and Erkoloc-pro have to be stored for app. 5 hours at a dry place without any pressure on it after having contact with alcohol in order to ensure that the alcohol can evaporate completely. Otherwise a bonding of the hard and soft layer is no longer guaranteed.

Sterilisation

A sterilisation with gaz and plasma (< 50 °C) is possible. As a result of the thermolability the materials are not autoclavable.

Mit freundlichen Grüßen

Erkodent Erich Kopp GmbH
Siemensstr.3
D-72285 Pfalzgrafenweiler
www.erkodent.com
info@erkodent.com
Tel.: +49 (0) 7445-8501-0
Fax: +49 (0) 7445-2092