
ERKODUR-A1

DIN EN ISO
13485: 2003

Erkoform-3/3d
Erkopress 300 Tp/-ci

0,6 mm, 170 °C
1,0 mm, 165 °C
2,0 mm, 155 °C

0,6 mm, 0:30 min
1,0 mm, 0:30 min
2,0 mm, 0:45 min

Erkoform-RVE
Erkopress 200 E

0,6 mm, 0:40 min
1,0 mm, 0:50 min
2,0 mm, 1:25 min

Lagerung • storage • stockage

ERKODENT Erich Kopp GmbH

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72285 Pfalzgrafenweiler, Germany

Tel. + 49 (0) 74 45/85 01-0

Fax + 49 (0) 74 45/20 92

info@erkodent.com • www.erkodent.com

ERKODENT®

DE ERKODUR-A1

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EN ERKODUR-A1

Chemical characteristics:

PETG: copolyester, PE, insulating foil: polyethylene
Harmless to health, BgVV (Federal Health Office)
listed thermoplasts with tested biocompatibility.
Insoluble in water, inactive, harmless to ground
water.

Disposal/Recycling:

Household waste (combustible, low in pollutants and
high in energy)/general recycling for plastic.

Field of use:

(see brochure thermoforming technique and materi-
al card) Stabilization and miniplast splints, tempo-
rary appliances, cosmetic splints.

Material properties:

Viscous hard foil. Extremely impact-resistant. Very
mouldable.
Bonds to acrylate (Resilit-S, 817 501 + 817 503).

Finishing:

(see brochure thermoforming technique) Finishing
set Quick 2 (110 877), high polish will be obtained
with polishes for plastic. The spacer foil will be taken
off only after thermoforming and finishing. Moistern
the area that has to be adjusted with acrylate at first
with little monomer.

Plastification:

(for machines without heating time control) Place
the foil in the devices in such a way that the
spacer/insulating foil is pointing towards the model.
Test softness of foil with an instrument. If permanent
impressions result, thermoform.

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