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Hints

Please observe mixing proportions.
Always keep Erkoton-HE/-W in a cool and dry place.
If colouring pastes are used, first mix material with colour pastes, then add hardener.
Mixing of Erkoton should be finished after 30 sec.
Warm material will accelerate setting time.
Erkoton combines very well to Sitran-N/F.
Lisko polishing discs, fine, white (413 678) are ideally suited for polishing the orthotic devices.
All numbers in brackets are order numbers.
Erkodent service phone + 49 (0) 74 45 / 85 01-45.

1. Erkoton can be mixed first without hardener paste to determine the required quantity. Excessive material that got in contact with hands cannot be put back into the pot.

4. When Erkoton and hardener paste have been roughly mixed, the mass is flattened in-between the palms until a uniform colour is achieved. Mixing of Erkoton should be finished within 30 sec.



7.-8. The ball or strip out of Erkoton will afterwards be formed directly on the patient's foot by smoothing or drawing. Firstly form the sole of foot (plantar). Then put the foot on the partial load plate.

10. The orthosis can now be formed under partial load (patient sits, foot stands on plate). Corrections on the toes must be held.



13.-14. Excessive parts can be removed with special scissors. After app. 12 min. the orthosis can be ground. Erkoton can be ground with pumicers, abrasive belts, multiple discs (Orthopedy) and cap grinders (Podology).


16. Use of colouring pastes

By colouring of Erkoton with colouring pastes the orthosis can be adjusted to skin colour or the pastes can serve as identifying feature for the patient (e.g. right and left).

18. Examples of different individually fabricated silicone orthoses with Erkoton-HE/-W:
a) interdigital toe wedge for Hallux valgus
b) pressure relief orthosis
c) toe support

2.-3. Press Erkoton flat in the hand and make an imprint into the material with the measuring spoon. For one level measuring spoon of Erkoton you require half a diameter of hardener paste, a strip of app. 2 cm. Mix Erkoton and hardener paste roughly with thumb or forefinger.

5.-6. A ball or strip (depending on the application) is formed out of the uniform mixed silicone material Erkoton in order to avoid air incorporations and folds during forming on the foot.

9. Now form the back of the foot (dorsal). By using the partial load plate both hands are free and can be used for correction and modelling of the orthosis.

11.-12. When producing a pressure relief orthosis, shortly before hardening of Erkoton a pressure relief can be worked in by pressing a finger into the material (no hole) to create a relief. After app. 8 min. Erkoton can be removed from the foot.

15. The surface is smoothed with Lisko polishing discs, fine, white.

17. The mixing of colouring pastes must be done before adding hardener paste. Please only use small quantities of colouring paste, as it is very intensive in colour.

The final hardening of Erkoton-HE/-W is reached after 24 hours. The Shore hardness of Erkoton-HE is app. 35 Shore A and of Erkoton-W app. 45 Shore A. Erkoton-HE/-W are medical products class I and correspond to the guidelines of Medical Devices.