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Hints

If SITRAN-N/F is coloured with colour pastes (red, yellow, brown) slightly more material should be mixed if the correction orthosis has to be drawn over again lateron.

The patient should be allowed to get used to the correction orthosis slowly. The correction orthosis should only be worn during the day.


LISKO polishing discs, fine, white are ideally suited for polishing the orthotic devices.

All numbers in brackets are article numbers.
ERKODENT service phone +49 74 45 / 85 01-45.


1.-2. Mix KNETON with hardener paste per instructions, this takes about 30 seconds. After about 3-4 min. remove the impression from the toe. The two side contours on the base of the impression must be extended by freshly mixed KNETON, so that the MOWAX model is provided with a base.


5.-6. The MOWAX model is placed in warm water at approx. 40 °C for 5 min. in order to soften it. After softening, the model is brought into the corrected position. Correct the model with a slight over-correction, in order to compensate the flexibility of the orthotic material.


9. The prepared toe model is then clamped in a small vice and a spacer drawn over. The spacer must sit tight without any air incorporations. The spacer enables a wearing of the orthosis the whole day and compensates swelling of the foot.


11. Cut tubular bandage double as long as the model and pull over isolated spacer and model base. The rest of the tubular bandage is put on and hangs backwards to lateron being easier turned and put over.

13.-14. Half of the mixed SITRAN-N/F is roughly applied to the tubular bandage so that the tubular bandage is moistened everywhere. The overhanging tubular bandage is turned around the own axis and pulled over the model moistened with SITRAN-N/F.

17.-18. Finally the insulated smoothing finger stall is slipped over. By this, the tubular bandage is moistened completely, the surface smoothed and the modelling of the mass is enabled. Apply SITRAN-N/F to correction and pressure relief spots and press thin in the area of transition.

21. After the final fitting the orthoses should be again coated with SITRAN-N/F in order to obtain a smooth surface (for this repeat steps 9+10, 17-21). Any possible air bubbles and the edges can be smoothed with LISKO polishing disc.

3.-4. Put a moulding support with the aid of a model holder in the KNETON negative and pour with MOWAX until the edge. After moulding wax MOWAX has solidified, the KNETON negative is opened with a wax knife and the MOWAX model removed.

7.-8. The model base is cut to size. The concave positions on the model base and unevenness (skin folds) can be smoothed out with liquid MOWAX. By removing the toe berry on the model, the correction can be raised.

10. The spacer is insulated with EKM powder. In the same operation the smoothing finger stall for smoothing the surface of the orthotic device is insulated on the inside and a small hole also pierced (so that air and excess material can escape).

12. All the required parts for example PeHa-Haft and smoothing finger stall must be prepared before the SITRAN-N/F is mixed as after mixture it remains onyl a 2 min. working time.

15.-16. The prepared PeHa-Haft is placed around the toe model tearing it slightly, the more PeHa-Haft is used, the more stable the orthosis will be. After that the remaining SITRAN-N/F should be applied mainly to the upper region, since it can then be more easily coated in a downwards direction afterwards.

19.-20. After app. 20 min. the smoothing finger stall can be removed. The orthoses is removed from the model and cut out roughly. The orthoses is then fitted on the patient in order to determine any possible pressure positions. Complete hardening takes place after app. 24 hours.

22. When fitting onto the patients, a little foot cream should be applied into the orthoses in order to facilitate slipping over. The correction orthoses can be cleaned with normal soap and water and can be treated with all commercial disinfectants.