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Hints

- Please observe the mixing proportions.
- Store Erkoton-20/30 at room temperature.
- Mixing of Erkoton should be finished after 30 sec.
- Warm material will accelerate setting time.
- Erkoton combines with Sitran-OF (413 056)
- Lisko polishing discs, fine, white (413 678) are ideally suited for polishing the orthotic devices.
- All numbers in brackets are order numbers.
- Erkodent service: info@erkodent.com

1. Take out the required quantity of Erkoton with the respective measuring spoons. **Please absolutely observe the markings of the pots and colours of the spoons.**


2. Take out Erkoton off the spoons and check quantity. **Mixing proportions 1 : 1.**

3. Combine Erkoton basic mass (brown) and catalyst (white) with each other and mix thoroughly in-between the finger tips.

4. After mixing Erkoton will be flattened in-between the palms to check the uniform colour. The mixing procedure should be finished within 30 sec.

5.-6. Depending on the application, form a ball or strip out of the uniformly mixed mass to avoid air incorporations and folds during forming. Put Erkoton on the desired position and form directly on the patient's foot by smoothing or drawing.


7. Firstly form the sole of foot (plantar). Then put the foot on the partial load plate in the podologique practice or on the measuring platform in the orthopaedic cabin.

8. Form the orthosis under partial load (patient sits, foot stands on plate). Corrections need to be held. By putting the foot on the plate or pedestal the hands can be used for correction and modelling. Now form the back of the foot (dorsal).



9. Correction needs to be held app. 3 min. after forming. When producing a pressure relief orthosis, shortly before hardening of Erkoton a pressure relief can be worked in by pressing a finger into the material (no hole) to create a relief. Pressure spots on the foot are thereby relieved.

10-11. In order to not destroy the detailed fitting, the Erkoton orthosis should not be removed before the expiration of the indicated time (app. 8 min.). It may be helpful to mark the form with a ballpen.

12. After app. 8 min. Erkoton can be removed from the foot.

13. Excessive parts can be removed with scissors, after app. 12 min. the orthosis can be ground. Erkoton can be ground with pumicers, abrasive belts, multiple discs (Orthopedy) and cap grinders (Podology).

14. The surface is smoothed with Lisko polishing discs, fine, white (413 678). The orthosis of Erkoton is now ready and can be worn directly on the foot.

15. Examples of different individually fabricated silicone orthoses with Erkoton 20/30:

- a) ball shell for Hallux valgus
- b) toe support
- c) interdigital toe wedge for Hallux valgus

16. The complete hardening of Erkoton 20/30 is reached after 24 h. The Shore hardness of Erkoton 20 is app. 20 Shore A and of Erkoton 30 app. 30 Shore A. Erkoton 20/30 are medical products class I and correspond to the guidelines of Medical Devices.