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Laboratory moulding material on divinyl polysiloxane based, addition curing, very high consistency, mouldable. Kneadable silicone filler to form sectional moulds for dental repairs in the laboratory, suitable when using cold and hot curing acrylic resins, etc. for isolation of artificial teeth, etc., boilable.

Read instructions for use

Keep away from sunlight

Keep dry

Recommended storage temperature

45 °C/113 °F

Recommended working temperature

25 °C/77 °F

15 °C/59 °F

## Dosage:

Mixture ratio: 1 : 1

One level measuring spoon (about 18 g) each catalyst (white) and base (yellow)

## Application:

Mix both components of **Aton Lab 80** about 30 sec. until the substance has a uniform colour.

## Hint:

Higher room temperatures will accelerate the setting.

## Technical Data:

- mixing can be done manually or by machine
- mixing time: **30 – 40 sec.**
- total working time: **1:30 min.**
- setting time: **4 min.**
- Shore A hardness: **approx. 80**

## Attention:

Take care that the catalyst (white) will not get dirty.

Do not mix with latex gloves.

Avoid stains on clothes.

***Please don't interchange spoons and lids!***

The user of our products is the only person responsible for the use and the application of our products. Our verbal and written advice has to be considered as advice without obligation also regarding possible protection rights of third parties. The user is not released from the obligation to check our products personally concerning suitability for the intended application. The perfect quality of our products is guaranteed according to our General Terms and Conditions. In any case and for all damage our liability is limited to the value of the goods supplied by us.

(manufacturer) Erkodent Erich Kopp GmbH

Siemensstraße 3 • 72285 Pfalzgrafenweiler, Germany • Tel. + 49 (0) 74 45/85 01-0

info@erkodent.com • www.erkodent.com • EN ISO 13485 • ISO 9001