Erkoplast PLA-R/-T/-W rosa/transparent/weiß

PLA-Griffe high/flat

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EN: Instructions for use

Intended use: Thermoforming sheets and accessories for the manufacture of apparatus for dental and orthodontic treatment.

Field of application: Erkoplast PLA-T/-W is thermoformed to fabricate intra-oral appliances such as individual impression trays, Erkoplast PLA-R is thermoformed to fabricate intra-oral appliances such as bite plates. PLA-handles to be attached to individual impression trays made from Erkoplast PLA.

Fabrication, see thermoforming technique brochure and videos, www.erkodent.com  

Contraindications: In case of already known allergic reaction to plastics a specific allergy to Erkoplast PLA-R/-T/-W (PLA) and PLA-handles should be excluded.

Precautionary measures: Erkoplast PLA-R/-T/-W and PLA-handles are only recommended for the above mentioned devices. Observe storage instructions.

Attention: The use is limited to the production of dental and orthodontic appliances. Use is the responsibility of the dentist. Use only on the recommendation or prescription of an authorised dentist. For single use only (thermoforming). Allergic reactions are unlikely, but cannot be excluded. Incorrect processing can lead to cracks, fractures and sharp edges. There is a risk of aspiration in the case of small breakages.

Erkoplast PLA-R/-T/-W (approx. values): density 1.24 g/cm³, modulus of elasticity 3500 MPa, elongation at break 5 %, water absorption 0.1 %, temp. resistance 55 °C

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