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**EN ISO • 13485:2016 • ISO 9001:2015**

## Instruction for use: **Erkosit**

**Labeling:**  consult accompanying documents  read instructions for use

 EC-directive 93/42/EEC

**Storage:**  keep away from sunlight  keep dry  95 °F  
35 °C  
5 °C  
41 °F recommended storage temperature

### Indications for use (general):

Self curing resin for the fabrication of dental appliances.

### Intended Use:

Erkosit is suitable for the production of temporary crowns and bridges.

### Kontraindikationen:

- Before use on patients with a history of allergic reactions to plastics it has to be clarified that there is no specific allergic reaction on Erkosit (polymethylmethacrylate).

### Warning:

- Use strictly limited for the fabrication of orthodontic and dental appliances.
- The use is subject to the responsibility of a therapist.
- For prescription use only.
- Allergic reactions are unlikely but possible.
- Hazardous substance according to EU CLP Regulation 1272/2008: GHS02, GHS07, GHS08

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- Erkosit is not recommended for other dental appliances than as described above under intended use.
- Hold safety data sheet for hazardous substances and notice (download, <https://www.erkodent.com>)
- Pay attention to the storage instructions.

### Chemical characteristics:

Self-curing polymethyl methacrylate based plastic, cadmium free, contains no tertiary amines.

**Disposal/Recycling:** General recycling for plastic if available, otherwise general waste.

### Material informations:

To the tooth color customizable liquid / powder system, repairable.

### Availabilities:

Erkosit set, each 100 g powder light / dark, 100 ml liquid 817 700 • 100 ml liquid 817 701  
100 g powder bright 817 702 • 100 g powder dark 817 703

### Technical data:

Vickers hardness 14 kp/mm<sup>2</sup> • Bending strength 75 N/mm<sup>2</sup> • Modulus of elasticity 2300/mm<sup>2</sup> • Impact strength 2.8 KJ/m<sup>2</sup>

**Working instruction:** Erkosit can be processed with the thermoforming foils Erkolen and Erkodur. Erkolen serves as a negative mold, which is removed after curing. Erkolen **does not bond** to Erkosit.

Erkodur in conjunction with Erkosit can be used to produce very stable temporary restorations. Erkodur **bonds to** Erkosit with the help of PV primer (222 970, 25 ml).

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### Product properties (at 21 °C):

Hardening system: self curing resin • Polymerization: approx. 12 min. • Processing time (mixing): approx. 2 min.

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Erkodent Erich Kopp GmbH • Siemensstraße 3 • 72285 Pfalzgrafenweiler, Germany