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 **ERKODENT**®

Instruction for use: **Resilit-S**

Labeling:  consult accompanying documents  read instructions for use

 EC-directive 93/42/EEC

Storage:  keep away from sunlight  keep dry  recommended storage temperature
5 °C / 41 °F 95 °F / 35 °C

Indications for use (general): Self curing resin for the fabrication of dental appliances.

Intended Use: Resilit-S is suitable for building up thermoformed splints, for the completion of model casting prostheses, for relining and prosthesis repairs.

Contraindication: Before use on patients with a history of allergic reactions to plastics it has to be clarified that there is no specific allergic reaction on Resilit-S (polymethylmethacrylate).

Warning:

- Use strictly limited for the fabrication of orthodontic and dental appliances.
- The use is subject to the responsibility of a therapist.
- Allergic reactions are unlikely but possible.
- Hazardous substance according to EU CLP Regulation 1272/2008: GHS02, GHS07, GHS08

Precautions:

- Resilit-S is not recommended for other dental appliances than as described above under intended use.
- Hold safety data sheet for hazardous substances and notice (Download, www.erkodent.com)
- Pay attention to the storage instructions.

Chemical characteristics: Self-curing polymethyl methacrylate based plastic, cadmium free, contains no tertiary amines.

 **Disposal/Recycling:** General recycling for plastic if available, otherwise general waste.

Material informations:

Clear or pink liquid / powder system, color-stable cold curing resin, high mechanical values.

Availabilities: 817 501 liquid 100 ml • 817 502 powder pink 100 g • 817 503 powder clear 100 g
817 601 liquid 500 ml • 817 602 powder pink 500 g • 817 603 powder clear 500 g

Technical data: Vickers hardness 16 kp/mm² • Bending strength 85 N/mm² • Modulus of elasticity 2500/mm²
water absorption 25.7 µg/mm³

Product properties (at 23 °C):

Hardening System: Self curing resin • Casting phase 2 - 3 minutes • Plastic phase 5 minutes

Working instruction: During liquid phase the RESILIT-S dough may be filled for example into a key, it no longer flows out of the key and can be modeled.

Moisten the area of the thermoformed splint to be built up with little monomer, then apply and model the Resilit-S dough.

After max. 8 minutes allow to polymerize in a pressure pot for 6 - 8 minutes at min. 2 bar and a water temperature of 40 - 50 ° C.

- Higher temperatures shorten, lower temperatures extend the aforementioned times.
- If needed, use insulating agent plaster against plastic.

 (manufacturer)

Erkodent Erich Kopp GmbH • Siemensstraße 3 • 72285 Pfalzgrafenweiler, Germany
Tel. + 49 (0) 74 45/85 01-0 • Fax + 49 (0) 74 45/85 01-15 • info@erkodent.com • www.erkodent.com

EN ISO • 13485:2016 • ISO 9001:2015

 **ERKODENT**®