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
**Ord. no.: 128 700**  
(With an additional  
fastening nut, 128 502)

Only one column?  
The upper and lower part are locked  
in place by a complex milled catch.  
At all times accurate identical initial  
and final positioning, set in stone.

Why only one column?  
Incomparably easy handling, free accessibility,  
best working control, fastest manufacturing of base  
and counterpart.

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
1. The enclosed guard tube protects the column and its toothing against impurities caused by plaster or resin. Remove the Sindos upper part, push the guard tube with the distance piece forward onto the column. The distance piece prevents that the guard tube seizes between plaster and column.

2. Pour the relining impression (generally the denture that has to be relined) with plaster. Put plaster onto the ground plate (lower part) and place the impression onto it. Pay attention that no plaster will get onto the toothing (or immediately remove it). Form the plaster that way that no undercuts will be generated, afterwards insulate (plaster against plaster, alginate based primer). Apply plaster for the counterpart onto the alignment.

3. + 4. Put the upper part on and seize it, pay attention to a gapless toothing. Fill counterpart with plaster through the retention opening and remove excessive material, allow to harden. Separate upper and lower part, the guard tube will stick either to the lower or upper part. Since the guard tube has some permissible variation, it can be loosened by compressing and easily taken off. Loosen the denture, remove the impression paste and roughen the inner surface of the denture.

5. Now push the protective guard with the distance piece backwards onto the lower column part. Spread the mixed resin over the plaster base.

6. Fix the denture in the counterpart and apply the mixed resin onto the roughened surface.

7. Put the upper part on, be careful to not bring any resin onto the toothing (or immediately remove it), and seize it, pay attention to a gapless toothing. Remove excessive resin.

For polymerization place Sindos into the pressure pot or polymerisation device (Liquisteam-E). Afterwards loosen the upper part and open the sindos. Knock the plaster out of the retention openings without damaging the device.

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Erkodent Erich Kopp GmbH • Siemensstraße 3 • 72285 Pfalzgrafenweiler • Germany