Erkoflex-95

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Erkoflex-95 (Ø 120 +125 mm und quadratisch 125 x 125 mm), transparent

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EN: Instructions for use

Erkoflex-95 (Ø 120 +125 mm and square 125 x 125 mm), transparent

The general purpose of Erkoflex-95 is the manufacture of appliances for dental and orthodontic treatment.

Application (material)	Intended use	Indication	Patient groups
Occlusal splints	CMD treatment	CMD	Adults, teenagers, children
Protective splints	Protection after local surgery	Local surgery (e.g. implantation)	Adults, teenagers

Contraindications: If an allergic reaction to plastics is already known, a specific allergy to Erkoflex-95 (EVA) should be ruled out.

Caution: Observe storage symbols.

User: Dental professionals  

Intended use: Erkoflex-95 is recommended as a medical device exclusively for the above-mentioned purposes.

Note: The use is subject to responsibility and is based on the recommendation/prescription of the dentist. Allergic reactions are unlikely, but cannot be ruled out.

Combination products: Erkoflex-95 can be processed in all devices of the Erkoform and Erkopress series from Erkodent. Erkoflex- 95 can be combined with Erkoflex Sticks-95.

Cleaning: Dental splints can be cleaned with Oxydens cleansing tablets. Also with curd soap, water and a toothbrush. Toothpaste is not suitable. Do not heat above 50 °C.

Please rinse the splint with clear water before inserting it. As the user, please inform the patient about the cleaning options.

Further information: Thermoforming technique brochure and videos, www.erkodent.com

Serious incidents in connection with Erkoflex-95 must be reported. The user must inform the competent authority and the manufacturer.

Performance characteristics (Circa values): Density 0.93 g/cm³, modulus of elasticity 70 MPa, elongation at break 450 %, water absorption 0.2 %, temp. resistance 62 °C

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