Erkoflex transparent und 1-farbig

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EN: Instructions for use

Erkoflex (round, Ø 120 +1 25 and square, 125 X 125 mm) transparent, 1-colored and colored as Erkoflex color

The general purpose of Erkoflex is the manufacture of appliances for dental and orthodontic treatment and for personal protective equipment.

Application (material)	Intended use	Indication	Patient groups
Radiation protection splint - Medical Device	Reduction of scattered radiation	malignant oral tumors	Adults, teenagers, children
Bracket transfer splint - MD	oral placement of brackets	misaligned teeth	Adults, teenagers
Multilayer protective splint - PPE	Bite and tooth protection	Risk of injury	Adults, teenagers, children
Positioner - MD	Positioning of teeth	Misaligned teeth	Adults, teenagers

Contraindications: If an allergic reaction to plastics is already known, a specific allergy to Erkoflex (EVA) should be ruled out.

Caution: Observe storage symbols.

User: Dental professionals

Intended use: Erkoflex is recommended exclusively for the above-mentioned purposes.

Note: The use is subject to responsibility and is based on the recommendation/prescription of the dentist. Allergic reactions are unlikely, but cannot be ruled out.

Combination products: Erkoflex can be processed in all devices of the Erkoform and Ekopress series from Erkodent. Erkoflex can be combined with Erkoflex Sticks-82.

Cleaning: Dental splints can be cleaned with Oxydens cleansing tablets. Also with curd soap, water and a toothbrush. Toothpaste is not suitable. Do not heat above 50 °C. Please rinse the splint with clear water before inserting it. As the user, please inform the patient about the cleaning options.

Further information: Thermoforming technique brochure and videos, www.erkodent.com

Serious incidents in connection with Erkflex must be reported. The user must inform the competent authority and the manufacturer.

Performance characteristics (Circa values): Density 0.95 g/cm³, modulus of elasticity 20 MPa, elongation at break 850 %, water absorption 0.5 %, temp. resistance 54 °C

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