

Product datasheet Erki set complete



1. Manufacturer information

Trade name:	Erki set complete
Intended use:	Nail correction
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

2. Intended use

Nail correction system out of plastic with different adjustable tensile forces, fits all nail forms.

3. Properties

Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	solid
Colour	-	transparent/ skin coloured
Odour	-	inodorous
Minimum durability	-	12 months

4. General information

Storage instructions:

Keep away from sunlight. Keep dry.

Recommended storage temperature: 5°C – 35°C

Recommended working temperature: 15°C – 25°C

Instructions for cleaning and maintenance:

Do not bring into contact with latex, disinfectants, foot creams or similar.

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.