

Product datasheet Erkoton 20/30



1. Manufacturer information

Trade name:	Erkoton 20/30
Intended use:	Fabrication of pressure relief orthoses
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

2. Intended use

Soft orthoses mass out of addition curing silicone for the fabrication of pressure relief orthoses directly on the patient foot.

3. Properties

Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	mouldable – type 0 – very high consistency
Colour	-	basic mass: brown
Odour	-	catalyst: grey
Density	ISO 1183	characteristic 1.22 g/cm ³
Hardness shore A	ISO 868	± 20/30
Minimum durability	-	18 months

Biological properties:

The material has been tested for biocompatibility according to DIN EN ISO 10993-1 and does not affect the patient's biological safety.

4. General information

Storage instructions:

Keep away from sunlight. Keep dry.

Recommended storage temperature: 5°C – 35°C

Recommended working temperature: 15°C – 25°C

Instructions for cleaning and maintenance:

Do not bring into contact with latex, disinfectants, foot creams or similar.

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.