

Product datasheet Mowax



1. Manufacturer information

Trade name:	Mowax
Intended use:	Modelling wax
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

2. Intended use

Modelling wax for the fabrication of correctable toe and forefoot models.

3. Composition

CAS-No.:	64748-51-4	63231-60-7
Designation:	Paraffin waxes	Hydrocarbon waxes, microcrystalline

4. Properties

Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	solid
Colour	-	red
Odour	-	characteristic
Density	ISO 1183	0.92 g/cm ³
Hardness shore A	-	± 44
Melting point	-	70°C

5. General information

Storage instructions:

Keep away from sunlight. Keep dry.

Recommended storage temperature: 5°C – 35°C

Recommended working temperature: 15°C – 25°C

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.