

Product datasheet Onycholit liquid



1. Manufacturer information

Trade name:	Onycholit liquid
Intended use:	Model nail plates
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

2. Intended use

For fixation of commercial wire clasps or for the fabrication of model nail plates for the Fraser clasp technique.

3. Properties

Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	liquid
Colour	-	colourless
Odour	-	characteristic, like: Ester
Density	ISO 1183	0.945 g/cm ³
Minimum durability	-	36 months

4. General information

Storage instructions:

Keep away from sunlight. Keep dry.

Recommended storage temperature: 5°C – 35°C

Recommended working temperature: 15°C – 25°C

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.