

## Product datasheet Onyclip



### 1. Manufacturer information

Trade name:	Onyclip
Intended use:	Nail correction
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

### 2. Intended use

Nail correction system out of high-grade steel with exactly adjustable tensile forces.

### 3. Properties

Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	coldhammered special steel, Material no.: 1.4310
Colour	-	white, beige
Strength	-	1500 – 1800 MPa
Odour	-	inodorous
Minimum durability	-	24 months

Elastic force of the Onyclip strips (1 cn ~ 1g)

Length [mm]	Spring travel [mm]	Elastic force [cn] (0,10)	Elastic force [cn] (0,15)
15	2	4	5
	4	7,2	9
20	2	3	4,2
	4	5	7,5

### 4. General information

Storage instructions super glue gel:

Keep away from sunlight. Keep cool and dry.

Storage life: app. 6 months after opening

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.