

Product datasheet super glue gel



1. Manufacturer information

Trade name:	Super glue gel
Intended use:	Super glue for nail correction
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

2. Intended use

Adhesive gel to fix Onyclip and Erki hooks.

3. Properties

Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	paste
Colour	-	colourless
Odour	-	pungent
Density	ISO 1183	1.09 g/cm ³
Minimum durability	-	12 months

4. General information

Storage instructions:

Keep away from sunlight. Keep dry.

Recommended storage temperature: 5°C – 35°C

Recommended working temperature: 15°C – 25°C

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.