

Product datasheet DIP paste



1. Manufacturer information

Trade name:	DIP paste
Intended use:	Paste for detection of potential sore spots
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

2. Intended use

Used to detect potential sore spots of new or existing dentures.

3. Properties

Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	paste
Colour	-	white
Odour	-	characteristic
Density	ISO 1183	> 0,9 g/cm ³
Minimum durability	-	24 months

Biological properties:

The material has been tested for biocompatibility according to DIN EN ISO 10993-1 and does not affect the patient's biological safety.

4. General information

Storage instructions:

Keep away from sunlight. Keep dry.
Recommended storage temperature: 5°C – 35°C
Recommended working temperature: 20°C – 22°C

Instructions for cleaning and maintenance:

Can be brushed off with warm water and soap.

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.