



## Erkoproess 240

Pressure forming unit for up to 6 models with base and for up to 9 models without base.

Erkoproess 240, 179 000



- The unit works with an overpressure of up to 4 bar. It thus meets the desire of many users to thermoform 3d print models with overpressure.
- The Erkoproess 240 enables the efficient production of correction splints/aligners or other series or larger numbers of splints.
- The thermoforming process is completely automated, the available plates are stored in the program and can be called up in the display with the respective parameters.
- Only thermoforming plates with a diameter of 240 mm can be processed.



Arrange models



Program machine



Forming in progress

