

## Product datasheet Resilit-S liquid



### 1. Manufacturer information

Trade name:	Resilit-S liquid
Intended use:	Repair and pouring resin
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

### 2. Intended use

Cold-curing acrylic resin to fabricate dentures, built up thermoformed splints, for the completion of model casting prostheses, for relining and prosthesis repairs.

### 3. Properties

#### Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	liquid
Colour	-	colourless
Odour	-	characteristic, like: ester
Density	ISO 1183	0,945 g/cm <sup>3</sup>
Flash point	-	10°C
Solvent content	-	0,0 %
Minimum durability	-	36 months

### 4. General information

#### Storage instructions:

Keep away from sunlight. Keep dry.

Recommended storage temperature: 5°C – 35°C

Recommended working temperature: 15°C – 25°C

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.