

## Product datasheet Wax bite rims



### 1. Manufacturer information

Trade name:	Wax bite rims
Intended use:	Fabrication of bite registrations
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

### 2. Intended use

Fabrication of bite registrations for full and partial prosthesis.

### 3. Properties

#### Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	solid
Colour	-	pink, green, orange, red
Odour	-	characteristic
Density	ISO 1183	0,92 g/cm <sup>3</sup> ± 34 (soft)
Hardness shore A	ISO 868	± 44 (medium) ± 55 (hard)
Minimum durability	-	36 months

#### Biological properties:

The material has been tested for biocompatibility according to DIN EN ISO 10993-1 and does not affect the patient's biological safety.

### 4. General information

#### Storage instructions:

Keep away from sunlight. Keep dry.  
Recommended storage temperature: 5°C – 35°C

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.